

# KAPA HyperCap SARS-CoV-2 panel



End to End viral surveillance workflow – from RNA to report

*May 2021*

*Markos Mihalatos, PhD – Sr. International Product Manager*



## Intended Use:

KAPA, HyperCap, HyperCapture, HyperPrep, HyperPlus and HyperDesign products are for Research Use Only. Not for use in diagnostic procedures.


Unless otherwise stated, all data generated by and on file with Roche.

HYPERCAP, HYPERCAPTURE, HYPERCHOICE, HYPEREXPLORE, HYPERDESIGN and KAPA are trademarks of Roche. All other product names and trademarks are property of their respective owners.

# SARS-CoV-2 spreads unprecedented challenges

First pandemic where the importance of Next Generation Sequencing (NGS) is broadly recognized

Roche



**genomeweb**

Business & Policy Technology Research Diagnostics Disease Areas Applied Markets Resources

### Biden Administration Will Provide \$1.6B to Improve COVID-19 Testing and Sequencing


Feb 17, 2021 | staff reporter

NEW YORK – The White House COVID-19 Response Team announced on Wednesday that the Biden administration will provide \$1.6 billion to expand and improve COVID-19 testing and genomic sequencing.

The US Department of Health and Human Services and the US Department of Defense will spend \$850 million to expand testing opportunities for schools kindergarten through eighth grade and underserved settings like homeless shelters.


The agency will also establish regional coordinating centers to identify laboratory testing capacity and match it to specific areas of need. The coordinating centers will partner with labs, including academic and commercial labs, to collect specimens, perform tests, and report results. Carole Johnson, the response team's supply coordinator, noted during a media briefing that testing can be difficult to implement in non-medical settings, so the coordinating centers will help bridge that gap.

B.1.1.7 and B.1.351 among SARS-CoV-2 variants \*



Adapted from GISAID

Spike mutations in B.1.351 (SA) and B.1.1.7 (UK) variants



- Ongoing **pandemic** has huge **health and economic cost**
- **Emerging variants** spread fast posing risks for vaccine efficiency
- The **need for fast and confident sequencing** is broadly acknowledged

## Genomic sequencing of SARS-CoV-2

A guide to implementation for maximum impact on public health

8 January 2021



Brussels, 19.1.2021  
COM(2021) 35 final

COMMUNICATION FROM THE COMMISSION TO THE EUROPEAN PARLIAMENT, THE EUROPEAN COUNCIL AND THE COUNCIL

A united front to beat COVID-19

# Sequencing is important in following viral spread and evolution

## Hybrid capture based target enrichment offers advantages for COVID-19 sequencing



Detect **phylogenetic relationships**  
(follow viral spread – surveillance and evolution)




Detect **nucleotide changes** in patient isolates  
(e.g. D614G, B.1.1.7, B.1.351)



Identify **potential new strains**



Detect **emerging resistance mutations** to  
antiviral drugs and vaccines



### Genomic sequencing of SARS-CoV-2

A guide to implementation for maximum impact on public health

8 January 2021



**“One advantage** of using a capture-based approach over a PCR amplicon-based approach (section 6.5.4) is that **capture-based** approaches can **tolerate sequence differences** from the probe sequences **of 10–20%**.

This is **higher than** the mismatch tolerated by PCR, where such a divergence from the primer sequences would result in a **high risk of amplicon failure.**”


<https://apps.who.int/iris/rest/bitstreams/1326052/retrieve>: 6.5.3 Targeted capture-based approaches

#### Sources:

James Hadfield, Colin Megill, Sidney M Bell, John Huddleston, Barney Potter, Charlton Callender, Pavel Sagulenko, Trevor Bedford, Richard A Neher, Nextstrain: real-time tracking of pathogen evolution, Bioinformatics, Volume 34, Issue 23, 01 December 2018, Pages 4121–4123, <https://doi.org/10.1093/bioinformatics/bty407>

# Viral surveillance with KAPA Sample Prep Tools and Exatype platform

*End-to-end workflow from RNA to analysis and meaningful reporting*



## KAPA RNA HyperCap Workflow for Sample Preparation with Exatype Platform analysis and reporting:

- **Single day & automatable** from RNA to Sequencer (with **1h hybridization**)
- **Single vendor offering** across the entire sample prep workflow
- **Long 120 bp probes** “resistant” to underlying variants
- **Secure on-cloud Exatype platform parallel analysis and reporting** by Hyrax Biosciences
- Purpose-built **user friendly analysis** with **5 simple steps** to **variant and lineage/clade reporting**

# Exatype Platform by Hyrax Biosciences

*Purpose built from viral sequence analysis experts for reliable SARS-CoV-2 surveillance*



High quality results

Confidently detect and report variants with the Exatype platform

Minimize time


Report your findings **faster by on-cloud parallel processing** using the intuitive interface of the Exatype platform

Scalability

Scale up with ease using Hyrax Biosciences' **streamlined and secure, cloud-based IT infrastructure**

Proven Expertise

Hyrax Biosciences' **proven expertise in viral sequence and reporting** offers **reliable SARS-CoV-2 surveillance**



Lineage/Clade Typing				
Sample	Pangolin (2.3.9 - 2021-04-23) <sup>1</sup>	Nextstrain clade (0.14.2) <sup>2</sup>	Nextstrain AA substitutions (0.14.2) <sup>2</sup>	Nextstrain AA deletions (0.14.2) <sup>2</sup>
Sample 1	B.1.1.7	✓ 20I/501YV1	ORF1a: T1001I, A1708D, I2230T ORF1b: P314L, T797I, K1383R S: N501Y, A570D, D614G, P681H, T716I, S982A, H1058N, D1118H ORF8: Q27*, R52I, Y73C N: M1X, D3L, R203K, G204R, S235F ORF7b: *44Q	S: Y144-

[DOWNLOAD AS CSV](#)  
<sup>1</sup><https://github.com/cov-lineages/pangolin>  
<sup>2</sup><https://clades.nextstrain.org/>



# Exatype Platform by Hyrax Biosciences

Purpose built from viral sequence analysis experts with intuitive user interface

Job details | Data format | Select data | Sample details | Review & submit

## Create a new job \*SARS-CoV-2

### Job details

Job Name  
SARS-CoV-2 surveillance ✓

Job Details (Optional- for your reference)  
Job Details (Optional- for your reference)

[NEXT >](#)



Job details | Data format | Select data | Sample details | Review & submit

## Create a new job \*SARS-CoV-2

### Data format

Are you submitting paired-end or single-ended data? Paired-end data requires two files per sample.

Illumina paired-end ✓  
 Illumina single-ended

[< PREV](#) [NEXT >](#)



Job details | Data format | Select data | Sample details | Review & submit

## Create a new job \*SARS-CoV-2

### Select data

Please select the data files that you would like to analyse. Only `gz`, `fastq` and `fq` files are accepted.

CHOOSE YOUR SEQUENCE DATA FILES

OR DRAG AND DROP YOUR SEQUENCE DATA FILES HERE

[< PREV](#) [NEXT >](#)

Job details | Data format | Select data | Sample details | Review & submit

## Create a new job \*SARS-CoV-2

### Sample details

#### Sample details | Illumina

Samples in this analysis [RESET SAMPLE IDS](#) [REMOVE ALL SAMPLES](#)

Sample ID	File name	Size	File status
Sample_100	Sample_100_R2.fastq.gz	0% 10.24 MB	Selected
	Sample_100_R1.fastq.gz	0% 9.60 MB	Selected

[< PREV](#) [NEXT >](#)



Job details | Data format | Select data | Sample details | Review & submit

## Create a new job \*SARS-CoV-2

### Review & submit

#### Job overview

Job name: SARS-CoV-2 surveillance Pathogen: SARS-CoV-2  
Number of samples: 1 Sequencing platform: Illumina  
Data structure: Paired-end Total data quantity: 19.84 MB

CLICK HERE TO UPLOAD 2 FILES AND RUN YOUR ANALYSIS 0/2 UPLOADED

Sample ID	File name	Size	File status
Sample_100	Sample_100_R2.fastq.gz	0% 10.24 MB	Selected
	Sample_100_R1.fastq.gz	0% 9.60 MB	Selected

[< PREV](#)



# Exatype Platform by Hyrax Biosciences

Purpose built from viral sequence analysis experts with intuitive user interface

Roche



My jobs

SHOW FILTERS

DOWNLOAD JOBS LOG

Page 1 of 1 30 items per page

Job name	Date	Info	Pathogen	Status	Platform	Analysis
SARS-CoV-2 surveillance	2021-05-09 06:14 pm	Show	SARS-CoV-2	Running 0%	Illumina	Pending

Analysis **starts automatically** after uploading



Lineage/Clade Typing

Sample	Pangolin (2.3.9 - 2021-04-23) <sup>1</sup>	Nextstrain clade (0.14.2) <sup>2</sup>	Nextstrain AA substitutions (0.14.2) <sup>2</sup>	Nextstrain AA deletions (0.14.2) <sup>2</sup>
UHR_100ng_MS2_Opg_SCV2_14_10000copy_3_DesignV3_singleplex_S21	B.1.1.7	✓ 20I/501Y.V1	ORF1a: T1001I, A1708D, I2230T ORF1b: P314L, T797I, K1383R S: N501Y, A570D, D614G, P681H, T716I, S982A, H1058N, D1118H ORF8: Q27*, R52I, Y73C N: M1X, D3L, R203K, G204R, S235F ORF7b: *44Q	S: Y144-

DOWNLOAD AS CSV


<sup>1</sup><https://github.com/cov-lineages/pangolin>  
<sup>2</sup><https://clades.nextstrain.org/>

Intuitive **reports generated automatically**




# Exatype Platform by Hyrax Biosciences

Purpose built from viral sequence analysis experts with intuitive user interface



Plethora of **downloadable reports**, from **sequencing metrics** through **variant frequency reports** and output files for **GISAID** and **Nextstrain**



### High-prevalence amino acid mutations

Observed at >=80% prevalence relative to the SARS-CoV-2 reference (MN908947.3)

Locus	Wild type	Locus position	Variant	Prevalence
ORF1a	T	1001	I	100%
ORF1a	A	1708	D	100%
ORF1a	I	2230	T	100%
ORF1b	P	314	L	100%
ORF1b	K	491	*	100%
ORF1b	T	797	I	100%
ORF1b	K	1383	R	100%
S	VY	143	V	100%
S	N	501	Y	100%
S	A	570	D	100%
S	D	614	G	100%
S	P	681	H	100%
S	T	716	I	100%
S	S	982	A	100%
S	H	1058	N	100%
S	D	1118	H	100%
ORF8	Q	27	*	100%
ORF8	R	52	I	100%
ORF8	Y	73	C	100%
N	D	3	L	100%
N	R	203	K	100%
N	G	204	R	100%
N	S	235	F	100%



### Low-prevalence amino acid mutations

Observed at <80% prevalence relative to the SARS-CoV-2 reference (MN908947.3)

⚠️ Low-prevalence mutations may be due to co-infection or within-host evolution, but also due to contamination. Proceed with caution.

Locus	Wild type	Locus position	Variant	Prevalence
ORF1a	L	3606	frameshift(TTtG)	22%
ORF1a	L	3606	frameshift(TTtG)	12%
ORF1a	Q	4392	K	11%
ORF1b	D	490	V	70%
ORF1b	I	2465	T	10%

### Downloads

All consensus sequences from SARS-CoV-2 samples

Combined

- Version3\_10k\_copies\_3\_consensus\_sequences.sars-cov-2.fasta

One-per-file

- Version3\_10k\_copies\_3\_aligned\_consensus\_sequences.zip

All mutations in SARS-CoV-2 samples present at >= 10% prevalence

- Version3\_10k\_copies\_3\_coding\_prevalences.zip
- Version3\_10k\_copies\_3\_noncoding\_prevalences.zip

Lineage/Clade files

- Version3\_10k\_copies\_3\_nextclade.csv
- Version3\_10k\_copies\_3\_nextclade.json
- Version3\_10k\_copies\_3\_nextclade.tree
- Version3\_10k\_copies\_3\_pangolin.csv

# KAPA Sample Prep Tools advance NGS COVID-19 research

*Benefits of KAPA RNA HyperPrep, KAPA HiFi DNA Polymerase and KAPA HyperCap Probes*



## KAPA RNA HyperPrep Kit including KAPA HiFi DNA Polymerase

Single-day automation-friendly library construction

KAPA HiFi DNA Polymerase

Robust and reliable performance

Higher success rates

inclusive of RNA enrichment

offers uniform and robust library amplification

across different sample types and input amounts

with lower input and degraded samples



## KAPA HyperCap Target Enrichment probes

Renowned probe design & content expertise

Manufactured with KAPA HiFi DNA Polymerase

High uniformity and low duplication

NGS probe pool QC

uncovers difficult genomic regions

enrich with higher fidelity and higher specificity

deliver high sequencing efficiency

provides consistent quality

# KAPA HyperCap SARS-CoV-2 Target Enrichment Design

*Efficient coverage of the viral genome with high performing probes*



## KAPA HyperCap SARS-CoV-2 Panel

- Covers **100% of the reference genome** (RefSeq) – target size ~30 Kb
- Covers >99.7% of another 183 genomic sequences from GenBank
- **Dense probe tiling**
- Probes targeting MS2 RNA for optional spike-in internal control
- **Based on hybrid capture** which **may advance** characterization, surveillance and viral evolution **research of divergent isolates**<sup>1</sup>

## Ordering details:


09436499001	KAPA HyperCap SARS-CoV-2 panel 12 rxn
09436502001	KAPA HyperCap SARS-CoV-2 panel 24 rxn
09436529001	KAPA HyperCap SARS-CoV-2 panel 48 rxn
09436537001	KAPA HyperCap SARS-CoV-2 panel 96 rxn

1. Genomic sequencing of SARS-CoV-2: a guide to implementation for maximum impact on public health. 8 Jan 2021 | *COVID-19: Laboratory and diagnosis*. World Health Organization.



# KAPA HyperCap SARS-CoV-2 Target Enrichment Panel

Research protocol combines KAPA RNA HyperPrep and HyperCap Probes benefits




# KAPA HyperCap SARS-CoV-2 Target Enrichment Panel

Research protocol combines KAPA RNA HyperPrep and HyperCap Probes benefits




rRNA depleted  
Or total RNA

KAPA RNA HyperPrep Kit





cDNA

KAPA RNA HyperPrep Kit



Pre-capture  
libraries

KAPA HyperCap v3.0 Workflow



## Minor protocol modifications

- Use 5  $\mu$ L of KAPA Universal Adapter in ligation (instead of 10  $\mu$ L)
- Omit first post ligation clean-up, perform only the 0.7X 2<sup>nd</sup> post ligation clean-up

# KAPA HyperCap SARS-CoV-2 Target Enrichment Panel


*Strong performance from low to high starting viral copies in the library preparation*

Roche

The KAPA HyperCap SARS-CoV-2 panel and target-enriched RNA-Seq workflow enables high-throughput sequencing of the SARS-CoV-2 genome

It was determined through testing the panel and workflow that:

- **About 97% 1x SARS-CoV-2 genome coverage** can be achieved **down to 1,000 viral copies** and 0.5 million NextSeq™ clusters (2 x 75 bp)
- **Genomic sequence (6% - 6.8%)** can be obtained from **as few as 10 viral input copies**
- SARS-CoV-2 coverage can further be increased by adding KAPA RiboErase for rRNA depletion prior to library preparation
- **Variants from six different SARS-CoV-2 isolates** were **detected** in the same sample **near their expected frequency**
- A **complex 10nt deletion** was detected **down to 1,000 total viral copies** (166 copies of the isolate with the deletion)



Data on file.

# KAPA HyperCap SARS-CoV-2 Target Enrichment Panel

## *Experimental overview and objectives*

**Demonstrate performance of the SARS-CoV-2 panel under the following conditions:**

- **From very low (10) to very high (1M) starting viral copies**
- **Pre-capture multiplexing of 8 or up to 16 sample libraries per hybridization**
- **One-hour compared to overnight hybridization**




Data on file.

# KAPA HyperCap SARS-CoV-2 Target Enrichment Panel




*Recovers viral reads even from just 10 starting viral copies*

- As few as **250k clusters** can cover the genome of just **10 viral copies** by **6% - 6.8%**
- Starting from **1000 viral copies**, **0.5M clusters** cover **92.4%** or **96.9%** of the genome in 20 ng or 100 ng UHR background



% bases  $\geq$  1x



Data on file.

20 ng

UHR background

100 ng




# KAPA HyperCap SARS-CoV-2 panel with pre-capture multiplexing

*Pre-capture multiplexing performance is concordant to single library capture*



Pre-capture multiplexing of 8 or up to 16 libraries (8-plex or 16-plex) compared to single captures (1-plex):


- Coverage metrics are similar across input viral copies for either 1-plex, 8-plex or up to 16-plex captures
- Robust capture performance offers convenience of pooling in the same hybridization



Data on file.

- When pre-capture multiplexing avoid mixing libraries with vastly different Ct values or viral loads
- **In the same capture high viral load libraries may compete for reads with the low viral load libraries**

### 8-plex Capture



Data on file.


20 ng UHRR, 100 ng UHRR

# KAPA HyperCap SARS-CoV-2 panel with 1-hour hybridization




Quick 1-hour hybridization offers similar performance to overnight in a single day workflow

## One hour hybridization offers good balance between workflow speed and performance at low viral concentrations



1 hour hyb samples subjected to 3 additional post-capture PCR cycles than overnight



20 ng UHRR background, 0.5 million NextSeq™ clusters (2 x 75 bp)

Data on file.


n=3 for each data point

# KAPA HyperCap SARS-CoV-2 Target Enrichment Panel



One hour hybridization provides broad coverage with low sequencing requirements

% bases uniquely covered at different levels vs 2 x 100 bp clusters from 10,000 viral copies



	12500	25000	50000	125000	250000	500000	1000000	4000000
1x	93,90%	98,03%	98,67%	99,00%	99,17%	99,23%	99,30%	99,40%
10x	2,30%	23,63%	66,17%	93,77%	97,37%	98,17%	98,33%	98,40%
20x	0,00%	0,43%	13,40%	63,57%	85,80%	92,63%	94,73%	95,80%
30x	0,00%	0,00%	1,00%	26,17%	56,60%	74,90%	81,53%	86,47%


## Robust performance with 1-hour hybridization:

- > 98% of the viral genome covered by 1x unique molecule depth with only 50,000 clusters of 2 x 100 bp.
- > 98% of the viral genome covered by 10x unique molecule depth with only 500,000 clusters of 2 x 100 bp.

Data on file. KAPA RNA HyperCap Workflow with 1 hour hybridization. Viral copies in 50 ng UHRR. Sequenced on a NovaSeq6000 System.

# KAPA HyperCap SARS-CoV-2 Target Enrichment Panel

*One hour hybridization provides broad coverage with low sequencing requirements*




## Robust performance with 1-hour hybridization:

- **23x mean unique coverage** with **only 125,000 clusters** of 2 x 100 bp.
- **39x mean unique coverage** with **only 500,000 clusters** of 2 x 100 bp.

Data on file. KAPA RNA HyperCap Workflow with 1 hour hybridization. Viral copies in 50 ng UHRR. Sequenced on a NovaSeq6000 System



# KAPA HyperCap SARS-CoV-2 minority variant detection

Concordant performance between expected and measured variant frequencies



1,000 viral copies (in 100 ng human RNA)

10,000 viral copies (in 100 ng human RNA)



Data on file.

Most complex variant – a 10nt deletion detected down to 1,000 and 10,000 copies in 100 ng human RNA

1 million NextSeq™ clusters (2 x 75 bp) – duplicates retained

# KAPA HyperCap SARS-CoV-2 Target Enrichment Panel


*Strong performance from low to high starting viral copies in the library preparation*

Roche

The KAPA HyperCap SARS-CoV-2 panel and target-enriched RNA-Seq workflow enables high-throughput sequencing of the SARS-CoV-2 genome


It was determined through testing the panel and workflow that:

- **About 97% 1x SARS-CoV-2 genome coverage** can be achieved **down to 1,000 viral copies** and 0.5 million NextSeq™ clusters (2 x 75 bp)
- **Genomic sequence (6% - 6.8%)** can be obtained from **as few as 10 viral input copies**
- SARS-CoV-2 coverage can further be increased by adding KAPA RiboErase for rRNA depletion prior to library preparation
- **Variants from six different SARS-CoV-2 isolates** were **detected** in the same sample **near their expected frequency**
- A **complex 10nt deletion** was detected **down to 1,000 total viral copies** (166 copies of the isolate with the deletion)



# KAPA Sample Prep Tools advance NGS COVID-19 research

*Breadth of additional offerings in the fight against the pandemic*





# Exatype Platform by Hyrax Biosciences

*Purpose built from viral sequence analysis experts for reliable SARS-CoV-2 surveillance*



High quality results

Achieve **uniform and broad genome coverage** with the **KAPA RNA HyperCap workflow** **Confidently detect and report variants** with the **Exatype platform**

Minimize time

**Save valuable time** and **maintain data quality** by **1-hour hybridization** with the KAPA HyperCap SARS-CoV-2 panel  
Report your findings **faster by on-cloud parallel processing** using the intuitive interface of the Exatype platform

Scalability


**Scale up with ease** using the **streamlined, automation friendly** KAPA HyperCap workflow and Hyrax Biosciences' **secure, cloud-based IT infrastructure**

Proven Technology

**Hybrid capture based methods may better detect relatively divergent SARS-CoV-2 sequences of new emerging variants<sup>1</sup>**

Proven Expertise

**Roche and Hyrax Biosciences** combine their **proven expertise in sample prep and reporting** to offer a **complete, end-to-end solution for SARS-CoV-2 surveillance**



# *Doing now what patients need next*


## **BACKUP SLIDES**

# KAPA HyperCap SARS-CoV-2 Target Enrichment Panel

*One hour hybridization provides broad coverage with low sequencing requirements*



% viral genome uniquely covered  $\geq 1x$  vs 2 x 100 bp clusters from 1,000 viral copies



**Robust performance with 1-hour hybridization from just 1,000 copies:**



- **More than 90% of the viral genome is uniquely covered by only 250,000 clusters of 2 x 100 bp.**

Data on file. KAPA RNA HyperCap Workflow with 1 hour hybridization. Viral copies in 50 ng UHRR. Sequenced on a NovaSeq6000 System

# KAPA HyperCap SARS-CoV-2 Target Enrichment Panel



*One hour hybridization provides broad coverage with low sequencing requirements*



## Robust performance with 1-hour hybridization with only 0.5M clusters of 2 x 100 bp:

- **> 99%** of the viral genome **uniquely covered** from **5,000 copies**
- **18x mean unique coverage** of the viral genome from **5,000 copies**

Data on file. KAPA RNA HyperCap Workflow with 1 hour hybridization. Viral copies in 50 ng UHRR. Sequenced on a NovaSeq6000 System

# KAPA HyperCap SARS-CoV-2 Target Enrichment Panel




*Great customer feedback for uniformity, specificity and genome coverage even at low virus loads*



Early Access Customer positive feedback:

- **Highly specific** (high on-target rates)
- **Highly sensitive** (reads from samples with Ct>32)
- **Broad Genome coverage**
- **Streamlined and user friendly** protocol



Data on file.

# KAPA HyperCap SARS-CoV-2 Target Enrichment Panel




*Great customer feedback for uniformity, specificity and genome coverage even at low virus loads*



Early Access Customer positive feedback:

- **Sequencing Uniformity**
- **Minority species analysis**
- **Cleaner data - No oligo-based amplification artifacts**
- **Streamlined protocol**

Normalized coverage of 1 (optimum)




Data on file.

# KAPA HyperCap SARS-CoV-2 Target Enrichment Panel

Offering key advantages over conventional PCR based approaches



## Workflow length comparisons



With RiboErase

Without RiboErase

Data on file.

## KAPA HyperCap SARS-CoV-2 with KAPA RNA HyperPrep:

- **Single tube** enrichment & automatable
- **Single day** from RNA to Sequencer
- **Long 120 bp probes** “resistant” to underlying variants
- **Single vendor offering** across the workflow
- **Superior uniformity** compared to amplicon based methods

## ARTIC based workflows:\*

- Legacy standard method
- Cannot call variants on either edge of the genome
- Split tubes (2) protocol
- New variants pose dropout risks, may need frequent updates
- Real-life TAT is 1.5 days to get to a sequencing library

\* <https://www.protocols.io/view/ncov-2019-sequencing-protocol-v3-locost-bh42j8ye>




# KAPA HyperCap SARS-CoV-2 compared to Competitor X panel

Strong performance against hybrid capture based competitor X at 250k clusters



Percent genome covered at 250k NextSeq™ clusters (2 x 75 bp)



## KAPA HyperCap SARS-CoV-2 panel

Provides better genome coverage with the same sequencing reads spent, than competitor X panel, at low or medium starting viral copy numbers


Data on file.

# KAPA HyperCap SARS-CoV-2 compared to Competitor X panel

*Strong performance against hybrid capture based competitor X at 0.5M clusters*



**Percent genome covered at 0.5M NextSeq™ clusters (2 x 75 bp)**



## **KAPA HyperCap SARS-CoV-2 panel**

Provides better genome coverage with the same sequencing reads spent, than competitor X panel, at low or medium starting viral copy numbers

Data on file.


# KAPA HyperCap SARS-CoV-2 panel with ribosomal RNA depletion



*RiboErase depletion improves coverage at low viral loads*

Depletion of rRNA with KAPA RiboErase provides better coverage at low viral concentrations

- **Broadening coverage** of the viral genome **from 92.3% to 98.8%**
- **Increasing ~3x** the percent of bases covered by  $\geq 10x$



Data on file.

100 ng UHRR background, 1 M NextSeq™ clusters (2 x 75 bp)